



## 吸头与移液器适配表

移液器品牌/型号/容量范围			普通吸头					滤芯吸头				
			0.5-10 $\mu$ l	10-200 $\mu$ l	10-300 $\mu$ l	100-1000 $\mu$ l	100-1250 $\mu$ l	0.5-10 $\mu$ l	10-100 $\mu$ l	10-200 $\mu$ l	100-500 $\mu$ l	100-1000 $\mu$ l
SOCOREX	CALIBRA/ACURA_Single channel	0.1-2 $\mu$ l	√					√				
		0.2-2 $\mu$ l	√					√				
		0.5-10 $\mu$ l	√					√				
		1-10 $\mu$ l	√					√				
		2-20 $\mu$ l		√								
		10-100 $\mu$ l		√					√			
		20-200 $\mu$ l		√	√					√		
	100-1000 $\mu$ l				√	√					√	
	CALIBRA/ACURA_Multi channel	1-10 $\mu$ l	√									
		5-50 $\mu$ l		√								
10-100 $\mu$ l			√									
20-200 $\mu$ l			√	√								
50-200 $\mu$ l			√	√			√				√	
VWR	SIGNATURE_Single channel	0.1-2 $\mu$ l	√					√				
		0.5-10 $\mu$ l	√					√				
		2-20 $\mu$ l		√					√			
		5-50 $\mu$ l		√					√			
		10-100 $\mu$ l		√						√		
		20-200 $\mu$ l		√	√					√		
		50-250 $\mu$ l		√	√				√			
	100-1000 $\mu$ l				√	√					√	
	SIGNATURE_Multi channel	10-100 $\mu$ l	√					√				
		5-50 $\mu$ l		√					√			
20-200 $\mu$ l			√	√					√			
50-300 $\mu$ l			√	√					√			
BIOHIT	PROLINE_Single channel	0.2-10 $\mu$ l	√					√	√			
		5-100 $\mu$ l		√					√			
		10-250 $\mu$ l		√	√					√		
		10-500 $\mu$ l				√	√				√	
		50-1000 $\mu$ l				√	√					√
	PROLINE_Multi channel	0.2-10 $\mu$ l	√					√				
		5-100 $\mu$ l		√					√			
		25-250 $\mu$ l		√	√					√		



## 吸头与移液器适配表

移液器品牌/型号/容量范围			普通吸头					滤芯吸头				
			0.5-10 $\mu$ l	10-200 $\mu$ l	10-300 $\mu$ l	100-1000 $\mu$ l	100-1250 $\mu$ l	0.5-10 $\mu$ l	10-100 $\mu$ l	10-200 $\mu$ l	100-500 $\mu$ l	100-1000 $\mu$ l
BIOHIT	E-LINE_Single channel	0.2-10 $\mu$ l	√					√				
		5-120 $\mu$ l		√					√			
		10-300 $\mu$ l		√	√					√		
		50-1000 $\mu$ l				√	√					√
	E-LINE_Multi channel	0.2-10 $\mu$ l	√					√				
		5-120 $\mu$ l		√					√			
		10-300 $\mu$ l		√	√							
	M-LINE_Single channel	0.5-10 $\mu$ l						√				
		2-20 $\mu$ l		√								
		10-100 $\mu$ l		√					√			
		20-200 $\mu$ l 100-1000 $\mu$ l		√	√		√	√				√
	M-LINE_Multi channel	0.5-10 $\mu$ l	√					√				
5-100 $\mu$ l			√					√				
30-300 $\mu$ l			√	√								
MATRIX	Impact/Impact2_Single channel	2-125 $\mu$ l		√						√		
		5-250 $\mu$ l								√		
		12-1250 $\mu$ l										
	Impact/Impact2_Multi channel	2-125 $\mu$ l		√							√	
		5-250 $\mu$ l 12-1250 $\mu$ l									√	
BRANDTECH	TRANSFERPETTE_Single channel	0.1-1 $\mu$ l	√					√				
		0.5-10 $\mu$ l	√					√				
		2-20 $\mu$ l		√								
		5-50 $\mu$ l		√					√			
		10-100 $\mu$ l		√					√			
		20-200 $\mu$ l		√	√					√		
		25-250 $\mu$ l 100-1000 $\mu$ l		√	√		√			√		
	TRANSFERPETTE_Multi channel	0.5-10 $\mu$ l	√					√				
		2-20 $\mu$ l		√				√				
		2.5-25 $\mu$ l		√				√				
		5-50 $\mu$ l		√								

## 吸头与移液器适配表

移液器品牌/型号/容量范围			普通吸头					滤芯吸头						
			0.5-10 $\mu$ l	10-200 $\mu$ l	10-300 $\mu$ l	100-1000 $\mu$ l	100-1250 $\mu$ l	0.5-10 $\mu$ l	10-100 $\mu$ l	10-200 $\mu$ l	100-500 $\mu$ l	100-1000 $\mu$ l		
BRANDTECH	TRANSFERPE TTE_Multi channel	10-100 $\mu$ l		√										
		20-200 $\mu$ l		√	√									
		30-300 $\mu$ l		√	√					√		√		
CAPP AERO	Pipettor_Single channel	0.2-10 $\mu$ l	√					√						
		5-10 $\mu$ l	√					√						
		2-20 $\mu$ l		√						√				
		10-100 $\mu$ l		√						√				
		25-200 $\mu$ l		√	√						√			
		50-300 $\mu$ l		√	√						√			
		100-1000 $\mu$ l				√	√					√		
	Pipettor_Multi channel	0.2-10 $\mu$ l	√					√						
		5-10 $\mu$ l	√					√						
		2-20 $\mu$ l		√						√				
		10-100 $\mu$ l		√						√				
		25-200 $\mu$ l		√	√						√			
		50-300 $\mu$ l		√	√						√			
		100-1000 $\mu$ l				√	√						√	
EPPENDORF	ELECTRONIC	0.5-10 $\mu$ l	√					√						
		5-100 $\mu$ l		√					√					
		20-200 $\mu$ l		√	√					√				
		50-1000 $\mu$ l				√	√					√		
	REFERENCE	0.1-1 $\mu$ l	√											
		2-20 $\mu$ l		√										
		10-100 $\mu$ l		√										
		20-200 $\mu$ l		√	√						√			
	RESEARCH_Si ngle channel	0.1-2.5 $\mu$ l	√											
		0.5-10 $\mu$ l	√											
		2-20 $\mu$ l		√										
		10-100 $\mu$ l		√					√					
		20- 200 $\mu$ l		√	√									
		100-1000 $\mu$ l				√	√						√	
		RESEARCH_M ulti channel	0.5-10 $\mu$ l	√					√					
			5-100 $\mu$ l		√					√				



## 吸头与移液器适配表

移液器品牌/型号/容量范围			普通吸头					滤芯吸头					
			0.5-10 $\mu$ l	10-200 $\mu$ l	10-300 $\mu$ l	100-1000 $\mu$ l	100-1250 $\mu$ l	0.5-10 $\mu$ l	10-100 $\mu$ l	10-200 $\mu$ l	100-500 $\mu$ l	100-1000 $\mu$ l	
EPENDORF	RESEARCH_Multi channel	50-100 $\mu$ l		√					√	√			
		20-300 $\mu$ l		√	√					√			
		30-300 $\mu$ l		√	√								
		50-1200 $\mu$ l				√	√						
THERMOLABSYSTEMS	FINNPIPETTE_Single channel	0.2-10 $\mu$ l	√					√					
		0.5-10 $\mu$ l	√					√					
		5-40 $\mu$ l		√					√				
		40-200 $\mu$ l		√	√					√			
		100-1000 $\mu$ l				√	√					√	
	200-1000 $\mu$ l				√	√					√		
	FINNPIPETTE_Multi channel	5-50 $\mu$ l		√					√				
		50-30 $\mu$ l		√	√					√			
GILSON	PIPETMAN_Single channel	0.1-2 $\mu$ l	√					√					
		0.2-10 $\mu$ l	√					√					
		2-20 $\mu$ l		√									
		10-100 $\mu$ l		√					√				
		20-200 $\mu$ l		√						√			
	100-1000 $\mu$ l		√	√							√		
PIPETMAN_Multi channel	20-200 $\mu$ l		√	√					√				
NICHIRYO/OFXO RD	Benchmate	0.5-10 $\mu$ l	√					√					
		2-20 $\mu$ l		√									
		10-50 $\mu$ l		√					√				
		10-100 $\mu$ l											
		20-200 $\mu$ l											
RAININ	EDP3_Single channel	0.5-10 $\mu$ l	√					√					
		2.5-25 $\mu$ l		√					√				
		10-100 $\mu$ l		√					√				
		25-250 $\mu$ l		√	√					√			
	100-1000 $\mu$ l				√	√					√		
	EDP3_Multi channel	0.5-10 $\mu$ l	√										
		10-250 $\mu$ l		√	√					√			
	LATCH MODE	0.1-2 $\mu$ l	√					√					
		0.2-10 $\mu$ l	√										
		2-20 $\mu$ l		√									
10-100 $\mu$ l			√					√					
20-200 $\mu$ l			√	√					√				
100-1000 $\mu$ l				√	√								