

ISO 21675-2019 配套标准品



百灵威供应最新全氟化合物标品定制，助力检测行业

ISO21675-PFAS：甲醇中 30 种全氟化合物

Analyte	CAS Number	Target Concentration
Perfluoro-n-butanoic acid	375-22-4	2 µg/mL
Perfluoro-n-pentanoic acid	2706-90-3	2 µg/mL
Perfluoro-n-hexanoic acid	307-24-4	2 µg/mL
Perfluoro-n-heptanoic acid	375-85-9	2 µg/mL
Perfluoro-n-octanoic acid	335-67-1	2 µg/mL
Perfluoro-n-nonanoic acid	375-95-1	2 µg/mL
Perfluoro-n-decanoic acid	335-76-2	2 µg/mL
Perfluoro-n-undecanoic acid	2058-94-8	2 µg/mL
Perfluoro-n-dodecanoic acid	307-55-1	2 µg/mL
Perfluoro-n-tridecanoic acid	72629-94-8	2 µg/mL
Perfluoro-n-tetradecanoic acid	376-06-7	2 µg/mL
Perfluorohexadecanoic acid	67905-19-5	2 µg/mL
Perfluorooctadecanoic acid (PFODA)	16517-11-6	2 µg/mL
Perfluorooctane sulfonamide	754-91-6	2 µg/mL
N-Methylperfluoro-1-octanesulfonamide	31506-32-8	2 µg/mL
Sulfuramid	4151-50-2	2 µg/mL
N-methylperfluoro-1-octanesulfonamidoacetic acid	2355-31-9	2 µg/mL
N-ethylperfluoro-1-octanesulfonamidoacetic acid	2991-50-6	2 µg/mL
2H-Perfluoro-2-decenoic Acid	70887-84-2	2 µg/mL
Perfluoro(2-methyl-3-oxahexanoic) acid	13252-13-6	2 µg/mL
Perfluorobutane-1-sulfonic acid	375-73-5	2 µg/mL
Perfluorohexane-1-sulfonic acid (Linear and Branched)	355-46-4	2 µg/mL
Perfluoroheptanesulfonic acid	375-92-8	2 µg/mL
Perfluorooctane-1-sulfonic acid	1763-23-1	2 µg/mL
Perfluorodecane-1-sulfonic acid	335-77-3	2 µg/mL
1H,1H,2H,2H-Perfluorooctane sulfonic acid	27619-97-2	2 µg/mL
1H,1H,2H,2H-Perfluorodecanesulfonic acid	39108-34-4	2 µg/mL
Sodium dodecafluoro-3H-4,8-dioxanonoate	N/A	2 µg/mL
Potassium 9-chlorohexadecafluoro-3-oxanone-1-sulfonate	73606-19-6	2 µg/mL
Bis[2-(perfluorooctyl)ethyl] phosphate	678-41-1	2 µg/mL



ISO 21675-2019 配套标准品



百灵威供应最新全氟化合物标品定制，助力检测行业

ISO 21675-LSS: ISO 21675 方法内标 (24 种)

ISO 21675-LSS	ISO 21675:2019 Labelled Stock Solution	1.2 mL
Perfluoro-n-(13C4)butanoic acid	MFFBA	100 ng/mL
Perfluoro-n-(13C5)pentanoic acid	M5PFPeA	100 ng/mL
Perfluoro-n-(1, 2, 3, 4, 6-13C5)hexanoic acid	M5PFHxA	100 ng/mL
Perfluoro-n-(1, 2, 3, 4-13C4)heptanoic acid	M4PFHpA	100 ng/mL
Perfluoro-n-(13C8)octanoic acid	M8PFOA	100 ng/mL
Perfluoro-n-(13C9)nonanoic acid	M9PFNA	100 ng/mL
Perfluoro-n-(1, 2, 3, 4, 5, 6-13C6)decanoic acid	M6PFDA	100 ng/mL
Perfluoro-n-(1, 2, 3, 4, 5, 6, 7-13C7)undecanoic acid	M7PFUdA	100 ng/mL
Perfluoro-n-(1, 2-13C2)dodecanoic acid	MPPDoA	100 ng/mL
Perfluoro-n-(1, 2-13C2)tetradecanoic acid	M2PFTeDA	100 ng/mL
Perfluoro-n-(1, 2-13C2)hexadecanoic acid	M2PFHxDA	100 ng/mL
Perfluoro-1-(13C8)octanesulfonamide	M8FOSA	100 ng/mL
N-methyl-d3-perfluoro-1-octanesulfonamide	d-N-MeFOSA	100 ng/mL
N-ethyl-d5-perfluoro-1-octanesulfonamide	d-N-EtFOSA	100 ng/mL
N-methyl-d3-perfluoro-1-octanesulfonamidoacetic acid	d3-N-MeFOSAA	100 ng/mL
N-ethyl-d5-perfluoro-1-octanesulfonamidoacetic acid	d5-N-EtFOSAA	100 ng/mL
2H-Perfluoro-(1, 2-13C2)-2-decenoic acid	MFOUEA	100 ng/mL
2, 3, 3, 3-Tetrafluoro-2-(1, 1, 2, 2, 3, 3, 3-heptafluoropropoxy) (13C3)propanoic acid	M3HFPO-DA	100 ng/mL
Sodium perfluoro-1-(2, 3, 4-13C3)butanesulfonate	M3PFBS	100 ng/mL
Sodium perfluoro-1-(1, 2, 3-13C3)hexanesulfonate	M3PFHxS	100 ng/mL
Sodium perfluoro-1-(13C8)octanesulfonate	M8PFOS	100 ng/mL
Sodium 1H, 1H, 2H, 2H-perfluoro(1, 2-13C2)octanesulfonate	M2-6:2FTS	100 ng/mL
Sodium 1H, 1H, 2H, 2H-perfluoro(1, 2-13C2)decanesulfonate	M2-8:2FTS	100 ng/mL
Sodium bis[1H, 1H, 2H, 2H-(1, 2-13C2)perfluorodecyl] phosphate	M4-8:2diPAP	100 ng/mL

