

Texas Red (λ_{ex} 596nm; λ_{em} 615nm) 检测通道, 荧光色: 红色

产品编号	产品名称	缩写	CAS	荧光性质
603231	Rhodamine B, 96% 罗丹明 B		81-88-9	λ_{ex} 553 nm; λ_{em} 627 nm in methanol
985419	Nile Red, 97.5% 尼罗红		7385-67-3	λ_{ex} 530 nm; λ_{em} 635 nm in methanol
279945	Propidium iodide, 95% 碘化丙啶	PI	25535-16-4	λ_{ex} 540 nm; λ_{em} 608 nm (bound to DNA)
1912458	Super Fluor 647, M.W. 900 - 1000		N/A	λ_{ex} 650 nm; λ_{em} 668 nm
137209	4-(4-Diethylaminostyryl)-N-methylpyridinium iodide, 95% 4-(4-二乙氨基苯乙烯基)-1-甲基碘化吡啶	4-Di-2-ASP	105802-46-8	λ_{ex} 485 nm; λ_{em} 607 nm in methanol
623461	Neurotrans Green C5, 95%		303727-66-4	λ_{em} 617 $\pm$ 4 nm in Methanol
2327737	Di6-ASP-BS, 95%		90133-77-0	λ_{em} 611 $\pm$ 4 nm in Methanol
208840	FM 1-43 荧光染料 FM 1-43 Dye, 95%		149838-22-2	λ_{em} 617 $\pm$ 4 nm in Methanol
220757	萘基荧光素 Naphthofluorescein, 95%	Naphthofluorescein	61419-02-1	λ_{ex} 594 nm; λ_{em} 663 nm in 0.1 M Tris pH 9.0
2327587	BDP 630/650 琥珀酰亚胺酯 BDP 630/650 NHS ester, 95%	BDP 630/650	380367-48-6	λ_{em} 635 $\pm$ 4 nm in Methanol
2327612	EverFluor TR-X, SE		217190-13-1	λ_{em} 616 $\pm$ 4 nm in Methanol
250555	4-(4-(二癸基氨基)苯乙烯)-N-甲基碘化吡啶 4-Di-10-ASP	4-Di-10-ASP	95378-73-7	λ_{em} 612 $\pm$ 4 nm in Methanol
594679	FM 2-10 荧光染料 FM 2-10 Dye, 95%		336185-20-7	λ_{ex} 505 $\pm$ 3 nm; λ_{em} 617 $\pm$ 4 nm in Methanol

620864	反-4-[4-(二甲氨基)苯乙烯基]-1-甲基 吡啶对甲苯磺酸 4-Di-1-ASP, TS		80969-52-4	$\lambda_{\text{em}} 605 \pm 4 \text{ nm}$
------------------------	--	--	------------	--